

Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit | CILP2 elisa kit

Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit

CATALOG #
AAA66916

REACTIVITY
Human


PRODUCT SYNONYMS

Cartilage Intermediate Layer Protein 2 (CILP2); Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit; AVSD2; CIRPIN; Atrioventricular Septal Defect 2; CILP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of CILP2. No significant cross-reactivity or interference between CILP2 and analogues was observed.
Sequence Length	1156
Samples	Tissue homogenates, cell lysates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.312-20ng/mL
Sensitivity	0.134ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level CILP2 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level CILP2 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CILP2 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of CILP2 in human tissue homogenates, cell lysates and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to CILP2. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to CILP2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain CILP2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of CILP2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	145699139
NCBI GeneID	148113
NCBI Accession #	NP_694953.2
NCBI GenBank Nucleotide #	NM_153221.2
UniProt Accession #	Q8IUL8; Q6NV88; Q8N4A6; Q8WV21
Molecular Weight	126,291 Da
NCBI Official Full Name	cartilage intermediate layer protein 2
NCBI Official Synonym Full Names	cartilage intermediate layer protein 2
NCBI Official Symbol	CILP2
NCBI Official Synonym Symbols	CLIP-2
NCBI Protein Information	cartilage intermediate layer protein 2
UniProt Protein Name	Cartilage intermediate layer protein 2
UniProt Gene Name	CILP2
UniProt Synonym Gene Names	CLIP-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cartilage-intermediate-layer-protein-2-cilp2/66916>