

Human Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) ELISA Kit | TIMP2 elisa kit

Human Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) ELISA Kit

CATALOG # AAA67105	REACTIVITY Human
------------------------------	----------------------------

PRODUCT SYNONYMS

Tissue Inhibitors Of Metalloproteinase 2 (TIMP2); Human Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) ELISA Kit; CSC-21K; TIMP Metallopeptidase Inhibitor 2; TIMP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	220
Samples	Serum, plasma, tissue homogenates and other biological fluids.
Assay Type	Sandwich
Detection Range	46.875-3000pg/mL
Sensitivity	14.5pg/mL
Intra Assay	CV<10%
Inter Assay	CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507511
NCBI GeneID	7077
NCBI Accession #	NP_003246.1
NCBI GenBank Nucleotide #	NM_003255.4
UniProt Accession #	P16035; Q16121; Q93006; Q9UDF7
Molecular Weight	24,399 Da
NCBI Official Full Name	metalloproteinase inhibitor 2
NCBI Official Synonym Full Names	TIMP metallopeptidase inhibitor 2
NCBI Official Symbol	TIMP2

NCBI Official Synonym Symbols	DDC8; CSC-21K
NCBI Protein Information	metalloproteinase inhibitor 2
UniProt Protein Name	Metalloproteinase inhibitor 2
UniProt Gene Name	TIMP2
UniProt Synonym Gene Names	TIMP-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tissue-inhibitors-of-metalloproteinase-2-timp2/67105>