

Human Decay Accelerating Factor (DAF) ELISA Kit | DAF elisa kit

Human Decay Accelerating Factor (DAF) ELISA Kit

CATALOG #
AAA67246

REACTIVITY
Human


PRODUCT SYNONYMS

Decay Accelerating Factor (DAF); Human Decay Accelerating Factor (DAF) ELISA Kit; CD55; CR; TC; Complement decay-accelerating factor; DAF elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	381
Samples	Serum, plasma and other biological fluids.
Assay Type	Sandwich
Detection Range	31.25-2000pg/mL
Sensitivity	13.7pg/mL
Intra Assay	CV<10%
Inter Assay	CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	10835143
NCBI GeneID	1604
NCBI Accession #	NP_000565.1
NCBI GenBank Nucleotide #	NM_000574.4
UniProt Accession #	P08174 ; P09679 ; P78361 ; Q14UF2 ; Q14UF3 ; Q14UF4 ; Q14UF5 ; Q14UF6 ; B1AP14 ; D3DT83 ; D3DT84 ; E7ER69
Molecular Weight	59,038 Da
NCBI Official Full Name	complement decay-accelerating factor isoform 1 preproprotein
NCBI Official Synonym Full Names	CD55 molecule (Cromer blood group)

NCBI Official Symbol	CD55
NCBI Official Synonym Symbols	CR; TC; DAF; CROM
NCBI Protein Information	complement decay-accelerating factor
UniProt Protein Name	Complement decay-accelerating factor
UniProt Gene Name	CD55
UniProt Synonym Gene Names	CR; DAF

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/decay-accelerating-factor-daf/67246>