

Horse Procollagen II N-Terminal Propeptide (PIINP) ELISA Kit | PIINP elisa kit

Horse Procollagen II N-Terminal Propeptide (PIINP) ELISA Kit

CATALOG #
AAA67279

REACTIVITY
Horse


PRODUCT SYNONYMS

Procollagen II N-Terminal Propeptide (PIINP); Horse Procollagen II N-Terminal Propeptide (PIINP) ELISA Kit; P2NP; N-Propeptide Of Type II Procollagen; Procollagen II Amino Terminal Propeptide; PIINP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Horse
Samples	Serum, plasma and other biological fluids.
Assay Type	Sandwich
Detection Range	3.125-200ng/mL
Sensitivity	1.34ng/mL
Intra Assay	CV<10%
Inter Assay	CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/horse-elisa-kits/procollagen-ii-n-terminal-propeptide-piinp/67279>