

Mouse Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit | BPI elisa kit

Mouse Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit

CATALOG # AAA67415	REACTIVITY Mouse
------------------------------	----------------------------


PRODUCT SYNONYMS

Bactericidal/Permeability Increasing Protein (BPI); Mouse Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit; CAP 57; BPI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	486
Samples	Serum, plasma and other biological fluids.
Assay Type	Sandwich
Detection Range	1.56-100pg/mL
Sensitivity	0.68pg/mL
Intra Assay	CV<10%
Inter Assay	CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	29244434
NCBI GeneID	329547
NCBI Accession #	NP_808518.1
NCBI GenBank Nucleotide #	NM_177850.3
UniProt Accession #	Q67E05 ; Q5I5I4 ; Q8BSF3 ; A2AC63
Molecular Weight	54,352 Da
NCBI Official Full Name	bactericidal permeability-increasing protein
NCBI Official Synonym Full Names	bactericidal permeability increasing protein

NCBI Official Symbol	Bpi
NCBI Official Synonym Symbols	Bpifd1; 9230105K17Rik
NCBI Protein Information	bactericidal permeability-increasing protein
UniProt Protein Name	Bactericidal permeability-increasing protein
UniProt Gene Name	Bpi
UniProt Synonym Gene Names	BPI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/bactericidal-permeability-increasing-protein-bpi/67415>