

Human Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) RTU ELISA Kit | ALPL rtu elisa kit

Human Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) ELISA Kit

CATALOG #
AAA69029

REACTIVITY
Human


PRODUCT SYNONYMS

Alkaline Phosphatase; Liver/Bone/Kidney (ALPL); Human Alkaline Phosphatase; Liver/Bone/Kidney (ALPL) ELISA Kit; BALP; BSAP; HOPS; AP-TNAP; TNSALP; Bone Alkaline Phosphatase; Tissue-Nonspecific Isozyme; ALPL rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	524
Assay Type	Sandwich
Samples	Serum, Plasma or Other Biological Fluids.
Detection Range	78.1-5000pg/mL
Sensitivity	36pg/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ALPL rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to ALPL. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to ALPL. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain ALPL, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of ALPL in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	116734717
NCBI GeneID	249
NCBI Accession #	NP_000469.3
NCBI GenBank Nucleotide #	NM_000478.5
UniProt Accession #	P05186; O75090; Q2TAI7; Q59EJ7; Q5BKZ5; Q5VTG5; Q6NZI8; Q8WU32; A1A4E7; B2RMP8; B7Z387; B7Z4Y6
Molecular Weight	51,045 Da
NCBI Official Full Name	alkaline phosphatase, tissue-nonspecific isozyme isoform 1 preproprotein
NCBI Official Synonym Full Names	alkaline phosphatase, liver/bone/kidney
NCBI Official Symbol	ALPL
NCBI Official Synonym Symbols	HOPS; TNAP; TNALP; APTNAP; TNSALP; AP-TNAP
NCBI Protein Information	alkaline phosphatase, tissue-nonspecific isozyme
UniProt Protein Name	Alkaline phosphatase, tissue-nonspecific isozyme
UniProt Gene Name	ALPL
UniProt Synonym Gene Names	AP-TNAP; TNSALP

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/alkaline-phosphatase-liver-bone-kidney-alpl/69029>