

Human Alkaline Phosphatase, Placental (ALPP) RTU ELISA Kit | ALPP rtu elisa kit

Human Alkaline Phosphatase, Placental (ALPP) ELISA Kit

CATALOG #
AAA69032

REACTIVITY
Human


PRODUCT SYNONYMS

Alkaline Phosphatase; Placental (ALPP); Human Alkaline Phosphatase; Placental (ALPP) ELISA Kit; PLAP1; ALP; ALPG; ALPPL; GCAP; Alkaline phosphatase Regan isozyme; Placental alkaline phosphatase 1; ALPP rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of ALPP. No significant cross-reactivity or interference between ALPP and analogues was observed.
Sequence Length	535
Samples	Serum, plasma, tissue homogenates or other biological fluids
Assay Type	Quantitative Sandwich
Detection Range	12.5-800pg/mL
Sensitivity	< 5.1pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level ALPP were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level ALPP were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ALPP rtu elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of ALPP in human serum, plasma, tissue homogenates or other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to ALPP. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to ALPP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain ALPP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of ALPP in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	94721246
NCBI GeneID	250
NCBI Accession #	NP_001623.3
NCBI GenBank Nucleotide #	NM_001632.4
UniProt Accession #	P05187; P05188; P06861; Q53S78; Q96DB7
Molecular Weight	57,954 Da
NCBI Official Full Name	alkaline phosphatase, placental type preproprotein
NCBI Official Synonym Full Names	alkaline phosphatase, placental
NCBI Official Symbol	ALPP
NCBI Official Synonym Symbols	ALP; PALP; PLAP; PLAP-1
NCBI Protein Information	alkaline phosphatase, placental type
UniProt Protein Name	Alkaline phosphatase, placental type
UniProt Gene Name	ALPP
UniProt Synonym Gene Names	PLAP; PLAP-1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/alkaline-phosphatase-placental-alpp/69032>