

## Mouse Carbohydrate Antigen 125 (CA125) RTU ELISA Kit | CA125 rtu elisa kit

Mouse Carbohydrate Antigen 125 (CA125) ELISA Kit

CATALOG #  
**AAA69470**

REACTIVITY  
**Mouse**


### PRODUCT SYNONYMS

Carbohydrate Antigen 125 (CA125); Mouse Carbohydrate Antigen 125 (CA125) ELISA Kit; MUC16; Cancer Antigen 125; Mucin 16; Cell Surface Associated; Ovarian Cancer-Related Tumor Marker CA125; Ovarian Carcinoma Antigen CA125; CA125 rtu elisa kit

### PRODUCT SPECIFICATIONS

|                           |                                                                                                                                                                                                                                                                                                                                                     |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                | Mouse                                                                                                                                                                                                                                                                                                                                               |
| Assay Type                | Sandwich                                                                                                                                                                                                                                                                                                                                            |
| Samples                   | Serum, Plasma, Tissue Homogenates or Other Biological Fluids.                                                                                                                                                                                                                                                                                       |
| Detection Range           | 0.312-20ng/mL                                                                                                                                                                                                                                                                                                                                       |
| Sensitivity               | 0.132ng/mL                                                                                                                                                                                                                                                                                                                                          |
| RTU ELISA Kit's Advantage | The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.                                                                                                                                                            |
| Preparation and Storage   | For unopened kits: All the reagents should be kept at 4 degree C upon receipt.<br>For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal. |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for CA125 rtu elisa kit

**Test Principle:** The microtiter plate provided in this kit has been pre-coated with an antibody specific to CA125. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to CA125. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain CA125, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm  $\pm$  10nm. The concentration of CA125 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

### NCBI AND UNIPROT INFORMATION

NCBI GeneID

[73732](#)

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://stage.aaabiotech.com/mouse-rtu-elisa-kits/carbohydrate-antigen-125-ca125/69470>