

Human Cartilage Intermediate Layer Protein 2 (CILP2) RTU ELISA Kit | CILP2 rtu elisa kit

Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit

CATALOG #
AAA69781

REACTIVITY
Human


PRODUCT SYNONYMS

Cartilage Intermediate Layer Protein 2 (CILP2); Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit; CILP2 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1156
Assay Type	Sandwich
Samples	Tissue Homogenates, Cell Lysates or Other Biological Fluids.
Detection Range	0.312-20ng/mL
Sensitivity	0.134ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CILP2 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to CILP2. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to CILP2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain CILP2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of CILP2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

145699139

NCBI GeneID	148113
NCBI Accession #	NP_694953.2
NCBI GenBank Nucleotide #	NM_153221.2
UniProt Accession #	Q8IUL8; Q6NV88; Q8N4A6; Q8WV21
Molecular Weight	126,291 Da
NCBI Official Full Name	cartilage intermediate layer protein 2
NCBI Official Synonym Full Names	cartilage intermediate layer protein 2
NCBI Official Symbol	CILP2
NCBI Official Synonym Symbols	CLIP-2
NCBI Protein Information	cartilage intermediate layer protein 2
UniProt Protein Name	Cartilage intermediate layer protein 2
UniProt Gene Name	CILP2
UniProt Synonym Gene Names	CILP-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/cartilage-intermediate-layer-protein-2-cilp2/69781>