

Human Eosinophil Peroxidase (EPX) RTU ELISA Kit | EPX rtu elisa kit

Human Eosinophil Peroxidase (EPX) ELISA Kit

CATALOG #
AAA70443

REACTIVITY
Human


PRODUCT SYNONYMS

Eosinophil Peroxidase (EPX); Human Eosinophil Peroxidase (EPX) ELISA Kit; EPO; EPP; EPX-PEN; EPER; EPX rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of EPX. No significant cross-reactivity or interference between EPX and analogues was observed.
Sequence Length	715
Assay Type	Quantitative Sandwich
Samples	Human serum, plasma, tissue homogenates or other biological fluids
Detection Range	1.25-80ng/mL
Sensitivity	< 0.51ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level EPX were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level EPX were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for EPX rtu elisa kit	Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to EPX. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to EPX. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain EPX, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of EPX in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of EPX in human serum, plasma, tissue homogenates or other biological fluids.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	4503595
NCBI GeneID	8288
NCBI Accession #	NP_000493.1
NCBI GenBank Nucleotide #	NM_000502.5
UniProt Accession #	P11678; Q4TVP3
Molecular Weight	81,040 Da
NCBI Official Full Name	eosinophil peroxidase preproprotein
NCBI Official Synonym Full Names	eosinophil peroxidase
NCBI Official Symbol	EPX
NCBI Official Synonym Symbols	EPO; EPP; EPXD; EPX-PEN
NCBI Protein Information	eosinophil peroxidase
UniProt Protein Name	Eosinophil peroxidase
UniProt Gene Name	EPX
UniProt Synonym Gene Names	EPER; EPO; EPP; EPO

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/eosinophil-peroxidase-epx/70443>