

Human Fatty Acid Desaturase 2 (FADS2) RTU ELISA Kit | FADS2 rtu elisa kit

Human Fatty Acid Desaturase 2 (FADS2) ELISA Kit

CATALOG #
AAA70531

REACTIVITY
Human


PRODUCT SYNONYMS

Fatty Acid Desaturase 2 (FADS2); Human Fatty Acid Desaturase 2 (FADS2) ELISA Kit; LLCDL2; FADSD6; D6D; TU13; DES6; Delta-6-Desaturase; Delta(6) fatty acid desaturase; FADS2 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of FADS2. No significant cross-reactivity or interference between FADS2 and analogues was observed.
Sequence Length	422
Assay Type	Quantitative Sandwich
Samples	Human serum, plasma, tissue homogenates or other biological fluids
Detection Range	78.1-5000pg/mL
Sensitivity	< 28pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level FADS2 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level FADS2 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FADS2 rtu elisa kit	Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to FADS2. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to FADS2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain FADS2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of FADS2 in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of FADS2 in human serum, plasma, tissue homogenates or other biological fluids.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	527498295
NCBI GeneID	9415
NCBI Accession #	NP_001268430.1
NCBI GenBank Nucleotide #	NM_001281501.1
UniProt Accession #	O95864; Q6MZQ7; Q96H07; Q96SV8; Q9H3G3; Q9Y3X4; A8K2M6; B7Z634
Molecular Weight	48,883 Da
NCBI Official Full Name	fatty acid desaturase 2 isoform 2
NCBI Official Synonym Full Names	fatty acid desaturase 2
NCBI Official Symbol	FADS2
NCBI Official Synonym Symbols	D6D; DES6; TU13; FADSD6; LLCDL2; SLL0262
NCBI Protein Information	fatty acid desaturase 2
UniProt Protein Name	Fatty acid desaturase 2
UniProt Gene Name	FADS2
UniProt Synonym Gene Names	D6D; Delta(6) desaturase; Delta-6 desaturase

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/fatty-acid-desaturase-2-fads2/70531>