

Mouse Glial Fibrillary Acidic Protein (GFAP) RTU ELISA Kit | GFAP rtu elisa kit

Mouse Glial Fibrillary Acidic Protein (GFAP) ELISA Kit

CATALOG #
AAA70879

REACTIVITY
Mouse


PRODUCT SYNONYMS

Glial Fibrillary Acidic Protein (GFAP); Mouse Glial Fibrillary Acidic Protein (GFAP) ELISA Kit; Intermediate Filament Protein; GFAP rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of GFAP. No significant cross-reactivity or interference between GFAP and analogues was observed.
Sequence Length	428
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids
Assay Type	Quantitative Sandwich
Detection Range	31.2-2000pg/mL
Sensitivity	< 10.0pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level GFAP were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level GFAP were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GFAP rtu elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of GFAP in mouse serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to GFAP. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to GFAP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain GFAP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of GFAP in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	196115327
NCBI GeneID	14580
NCBI Accession #	NP_001124492.1
NCBI GenBank Nucleotide #	NM_001131020.1
UniProt Accession #	P03995; Q09J71; Q3USS4; Q496R4; Q496S3; Q7TQ30; Q80VX6; Q925K2; Q925K3; A1E2H7; A2AH87; B2RTI7
Molecular Weight	49,365 Da
NCBI Official Full Name	glial fibrillary acidic protein isoform 1
NCBI Official Synonym Full Names	glial fibrillary acidic protein
NCBI Official Symbol	Gfap
NCBI Official Synonym Symbols	A1836096
NCBI Protein Information	glial fibrillary acidic protein
UniProt Protein Name	Glial fibrillary acidic protein
UniProt Gene Name	Gfap
UniProt Synonym Gene Names	GFAP

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-rtu-elisa-kits/glial-fibrillary-acidic-protein-gfap/70879>