

Human Glutamate Receptor, Ionotropic, AMPA 1 (GRIA1) RTU ELISA Kit | GRIA1 rtu elisa kit

Human Glutamate Receptor, Ionotropic, AMPA 1 (GRIA1) ELISA Kit

CATALOG # AAA71008	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

Glutamate Receptor; Ionotropic; AMPA 1 (GRIA1); Human Glutamate Receptor; AMPA 1 (GRIA1) ELISA Kit; GLUH1; GLUR1; GLURA; HBGR1; GluR-K1; AMPA-selective glutamate receptor 1; GRIA1 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of GRIA1. No significant cross-reactivity or interference between GRIA1 and analogues was observed.
Sequence Length	906
Samples	Human Tissue Homogenates Or Other Biological Fluids
Assay Type	Quantitative Sandwich
Detection Range	0.156-10ng/mL
Sensitivity	0.056ng/mL.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level GRIA1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level GRIA1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GRIA1 rtu elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of GRIA1 in human tissue homogenates or other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to GRIA1. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to GRIA1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain GRIA1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of GRIA1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	167001419
NCBI GeneID	2890
NCBI Accession #	NP_000818.2
NCBI GenBank Nucleotide #	NM_000827.3
UniProt Accession #	P42261; Q2NKM6; B7Z2S0; B7Z2W8; B7Z3F6; B7Z9G9; D3DQI4; E7ESV8
Molecular Weight	102,864 Da
NCBI Official Full Name	glutamate receptor 1 isoform 1
NCBI Official Synonym Full Names	glutamate ionotropic receptor AMPA type subunit 1
NCBI Official Symbol	GRIA1
NCBI Official Synonym Symbols	GLUH1; GLUR1; GLURA; GluA1; HBGR1
NCBI Protein Information	glutamate receptor 1
UniProt Protein Name	Glutamate receptor 1
UniProt Gene Name	GRIA1
UniProt Synonym Gene Names	GLUH1; GLUR1; GluR-1; GluA1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.