

Human Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) RTU ELISA Kit | GRIA2 rtu elisa kit

Human Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) ELISA Kit

CATALOG # AAA71010	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

Glutamate Receptor; Ionotropic; AMPA 2 (GRIA2); Human Glutamate Receptor; AMPA 2 (GRIA2) ELISA Kit; GLUR2; GLURB; HBGR2; GluR-K2; AMPA-selective glutamate receptor 2; GRIA2 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	883
Assay Type	Sandwich
Samples	Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids.
Detection Range	0.312-20ng/mL
Sensitivity	0.115ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GRIA2 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to GRIA2. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to GRIA2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain GRIA2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of GRIA2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	74275354
NCBI GeneID	2891
NCBI Accession #	NP_000817.2
NCBI GenBank Nucleotide #	NM_000826.3
UniProt Accession #	P42262; Q96FP6; A8MT92; I6L997
Molecular Weight	93,775 Da
NCBI Official Full Name	glutamate receptor 2 isoform 1
NCBI Official Synonym Full Names	glutamate ionotropic receptor AMPA type subunit 2
NCBI Official Symbol	GRIA2
NCBI Official Synonym Symbols	GLUR2; GLURB; GluA2; HBGR2; GluR-K2
NCBI Protein Information	glutamate receptor 2
UniProt Protein Name	Glutamate receptor 2
UniProt Gene Name	GRIA2
UniProt Synonym Gene Names	GLUR2; GluR-2; GluA2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/glutamate-receptor-ionotropic-ampa-2-gria2/71010>