

Nitric Oxide Non-Enzymatic Assay Kit

Nitric Oxide Non-Enzymatic Assay Kit

CATALOG #
AAA71820


PRODUCT SYNONYMS

Nitric Oxide Non-Enzymatic; Nitric Oxide Non-Enzymatic Assay Kit; Non-Enzymatic Nitric Oxide Assay Kit; Nitric Oxide Non-Enzymatic assay kit

PRODUCT SPECIFICATIONS

Assay Type	Sandwich
Preparation and Storage	1. The Cadmium Beads should be stored in a dedicated dry environment. 2. Store this kit and its components at 4 degree C until use. 3. Do not freeze..

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Nitric Oxide Non-Enzymatic assay kit	Principle of the Assay: Human Total Prostate Specific Antigen (tPSA) ELISA Kit employs a two-site sandwich ELISA to quantitate tPSA in samples. An antibody specific for Human tPSA has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any tPSA present is bound by the immobilized antibody. After removing any unbound substances, HRP-Conjugated tPSA detection antibody is added to the wells. Following a wash to remove any unbound HRP reagent, a Chromogen solution is added to the wells and color develops in proportion to the amount of tPSA bound in the initial step. The color development is stopped and the intensity of the color is measured.
References	1. Schmidt, H.H., et. al., (1995) Biochemica 2:22-23

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/assay-kits/nitric-oxide-non-enzymatic/71820>