

Goat anti Bovine IgG (FITC) Secondary Antibody | IgG secondary antibody

Goat anti Bovine IgG (FITC)

CATALOG # AAA75476	HOST Goat	REACTIVITY Bovine IgG	APPLICATIONS Immunofluorescence, Flow Cytometry, Functional Assay, ELISA
------------------------------	---------------------	---------------------------------	--


PRODUCT SYNONYMS

Goat anti Bovine IgG (FITC); Goat anti Bovine Immunoglobulin G secondary (FITC); IgG secondary antibody

PRODUCT SPECIFICATIONS

Host	Goat
Reactivity	Bovine IgG
Specificity	Bovine IgG
Purity/Purification	Goat anti Bovine IgG (FITC) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.
Form/Format	Lyophilized from 0.01M Na3PO4, pH 7.4, with 0.14M NaCl, 10 mg/ml BSA, and 0.01% NaN3. Immunoglobulin and protease free.
Applicable Applications for IgG secondary antibody	IF (Immunofluorescence), FCM/FACS (Flow Cytometry), ELISA
Biohazard Information	This product contains sodium azide as preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.
Tag/Conjugate	FITC
Immunogen	Goat anti-bovine IgG (FITC) was raised in goat using bovine IgG F(ab')2 fragment as the immunogen.
Preparation and Storage	Store at 4 degree C until reconstitution. Following reconstitution aliquot and freeze at -20 degree C for long term storage. Avoid repeated freeze/thaw cycles.

IMAGES



DB (Dot Blot)

Dot Blot of Goat anti-Bovine IgG Fab2 Antibody Fluorescein Conjugated. Antigen: Bovine IgG.

Load: Lane 1 - 100 ng Lane 2 - 33.3 ng Lane 3 - 11.1 ng Lane 4 - 3.70 ng Lane 5 - 1.23 ng.

Primary antibody: n/a. Secondary antibody: Goat anti-Bovine IgG Fab2 Antibody Fluorescein

Conjugated at 1:1,000 for 1 HR at RT. Block: MB-070 for 1 HR at RT.

ADDITIONAL INFORMATION

Related Product Information for IgG secondary antibody	Goat anti Bovine IgG secondary antibody (FITC)
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/goat-bovine-ovine-secondary-antibody/goat-anti-bovine-igg-fitc/75476>