

Goat anti Human IgG (H + L) (FITC) Secondary Antibody | IgG secondary antibody

Goat anti Human IgG (H + L) (FITC)

CATALOG #	HOST	REACTIVITY	APPLICATIONS
AAA75619	Goat	Assay by immunoelcerophoresis resulted in a single precipitin arc against anti-Fluorescein, anti-goat serum, Human serum.	Immunofluorescence, Flow Cytometry, Functional Assay, ELISA


PRODUCT SYNONYMS

Goat anti Human IgG (H + L) (FITC); Goat anti Human Immunoglobulin G (H + L) secondary (FITC); IgG secondary antibody

PRODUCT SPECIFICATIONS

Host	Goat
Reactivity	Assay by immunoelcerophoresis resulted in a single precipitin arc against anti-Fluorescein, anti-goat serum, Human serum.
Specificity	Human IgG (H + L)
Purity/Purification	Goat anti Human IgG (H + L) (FITC) was purified by affinity chromatography.
Form/Format	Liquid in 0.01M Sodium Ohosphate, 0.15M Sodium Chloride, pH 7.2, with 10mg/ml BSA and 0.01% Thimerosal. Reconstitute with 1.0ml deionized water.
Applicable Applications for IgG secondary antibody	IF (Immunofluorescence), FCM/FACS (Flow Cytometry), ELISA
Tag/Conjugate	FITC
Immunogen	Goat anti-human IgG (H+L) (FITC) was raised in goat using human IgG whole molecule as the immunogen.
Biological Significance	Immunoglobulin G (IgG) are antibody molecules. Each IgG is composed of four peptide chains - two heavy chains ? and two light chains. Each IgG has two antigen binding sites. Other Immunoglobulins may be described in terms of polymers with the IgG structure considered the monomer.
Preparation and Storage	Ships ambient or refrigerated. Upon receipt store at 4°C prior to restoration. Once reconstituted aliquot and freeze at -20°C. Avoid repeat freeze/thaw cycles.

IMAGES



WB (Western Blot)

ADDITIONAL INFORMATION

Related Product Information for IgG secondary antibody	Goat anti Human IgG (H + L) secondary antibody (FITC)
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/goat-goat-human-secondary-antibody/goat-anti-human-igg-h-l-fitc/75619>