

Agarose Gel DNA Extraction and Purification Extraction Kit

Agarose Gel DNA Extraction and Purification

CATALOG #
AAA76700

PRODUCT SYNONYMS


Agarose Gel DNA Extraction and Purification; Agarose Gel DNA Extraction and Purification extraction kit

PRODUCT SPECIFICATIONS

Preparation and Storage

Mixed storage temperatures. Magnetic beads stored at 2-8 degree C. Other reagents stored at room temperature. Avoid repeated freeze-thaw cycles.
Shipping: Refrigerated (polar packs)

IMAGES



Application Data

ADDITIONAL INFORMATION

Related Product Information for Agarose Gel DNA Extraction and Purification extraction kit

This kit provides a simple, rapid and efficient method for the recovery and purification of DNA directly from Agarose gel (100 bp to 50 kb) with typical recovery efficiency up to 85%. The resulting product can be directly used for sequencing, restriction digestion, or PCR and other downstream experiments. In addition, the kit can be used to concentrate DNA. The kit will work with a 96 well round bottom plates if a special magnetic frame is used. The kit can also be used with a variety of automatic nucleic acid extraction instruments and workstation.

Product Categories/Family for Agarose Gel DNA Extraction and Purification extraction kit

[Kits](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/extraction-kits/agarose-gel-dna-extraction-and-purification/76700>