

Si-Mag Magnet Rack | Si-Ma rack

Si-Mag Magnet

CATALOG #
AAA76702

PRODUCT SYNONYMS

Si-Mag Magnet; Si-Ma rack

IMAGES



Application Data



Application Data

ADDITIONAL INFORMATION

Related Product Information for Si-Ma rack

Multi-function strong magnet rack to be used with Si-Mag bead Kits
Ergonomic racks for efficient separation of nucleic acids, antibodies, proteins and other biomolecules using magnetic beads and standard plastic tubes of three different sizes
Strong magnets embedded into the removable block for fast and efficient sequestration and immobilization of magnetic beads with material bound to them
Magnet block can be removed from the rack for mixing samples in tubes inserted into the rack
In horizontal position holds up to 12 standard laboratory 1.5 - 2 mL tubes (top image) and in vertical position holds one 15 mL and one 50 mL tube (bottom image)
Unique design of collar tube clamps: tubes with samples are firmly held in place during mixing samples by flipping rack upside down
Contents of tubes and their labels can be easily seen from different angle views
Durable, compact and comfortable to handle
Easy to clean and disinfect

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/si-ma-rack/si-mag-magnet/76702>