

Chicken Campylobacter Jejuni Polyclonal Antibody

Anti-Campylobacter Jejuni Antibody-Affinity Purified

CATALOG #	HOST	REACTIVITY	APPLICATIONS
AAA76869	Chicken	Campylobacter jejuni	Western Blot, Immunohistochemistry, Immunofluorescence, ELISA

PRODUCT SYNONYMS

Campylobacter Jejuni; Polyclonal Antibody; Anti-Campylobacter Jejuni Antibody-Affinity Purified; anti-Campylobacter Jejuni Antibody (HipO); anti-Campylobacter Jejuni antibody


PRODUCT SPECIFICATIONS

Host	Chicken
Reactivity	Campylobacter jejuni
Clonality	Polyclonal
Isotype	IgY
Specificity	This antibody reacts with Campylobacter jejuni hippurate hydrolase as determined by ELISA and it can detect Campylobacter jejuni cells with Immunofluorescence staining.
Purity/Purification	Antigen Affinity Purified
Form/Format	Liquid; Phosphate buffered saline, pH7.2; 0.09% sodium azide
Applicable Applications for anti-Campylobacter Jejuni antibody	WB (Western Blot), IHC (Immunohistochemistry), IF (Immunofluorescence), ELISA
Type	Unconjugated AP
Sample Type	Control
Procedure	Chickens were immunized with highly purified Campylobacter jejuni hippurate hydrolase (HipO) and the eggs were collected and delipidated. Antibody was immunoaffinity purified off a HipO containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.
Preparation and Storage	Store at 2-8 degree C Shelf life: 24 months from shipping date

IMAGES



IF (Immunofluorescence)
Affinity Purified Chicken anti-Campylobacter Jejuni (202)-Reacts with Hippurate hydrolase.



IF (Immunofluorescence)
Affinity Purified Chicken anti-Campylobacter Jejuni (202)-Reacts with Hippurate hydrolase.

ADDITIONAL INFORMATION

Product Categories/Family for anti-Campylobacter Jejuni antibody	Primary Antibodies; Infectious Agents-Human and Veterinary; Campylobacter jejuni Antibodies
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/polyclonal/campylobacter-jejuni-antibody-western-blot-immunohistochemistry-immunofluorescence-elisa/76869>