

Chicken Campylobacter Jejuni Polyclonal Antibody

Anti-Campylobacter Jejuni Antibody-Affinity Purified

CATALOG # AAA76871	HOST Chicken	REACTIVITY Campylobacter jejuni	APPLICATIONS Western Blot, Immunohistochemistry, Immunofluorescence, ELISA
------------------------------	------------------------	---	--


PRODUCT SYNONYMS

Campylobacter Jejuni; Polyclonal Antibody; Anti-Campylobacter Jejuni Antibody-Affinity Purified; anti-Campylobacter Jejuni Antibody (CadF); anti-Campylobacter Jejuni antibody

PRODUCT SPECIFICATIONS

Host	Chicken
Reactivity	Campylobacter jejuni
Clonality	Polyclonal
Isotype	IgY
Specificity	This antibody reacts with Campylobacter jejuni fibronectin binding outer membrane protein as determined by ELISA and it can detect Campylobacter jejuni cells with Immunofluorescence staining.
Purity/Purification	Antigen Affinity Purified
Form/Format	Liquid; Phosphate buffered saline, pH7.2; 0.09% sodium azide
Applicable Applications for anti-Campylobacter Jejuni antibody	WB (Western Blot), IHC (Immunohistochemistry), IF (Immunofluorescence), ELISA
Type	Unconjugated AP
Sample Type	Control
Procedure	Chickens were immunized with highly purified Campylobacter jejuni Fibronectin binding outer membrane protein (CadF) and the eggs were collected and delipidated. Antibody was immunoaffinity purified off CadF containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.
Preparation and Storage	Store at 2-8 degree C Shelf life: 24 months from shipping date

IMAGES



IF (Immunofluorescence)
Affinity Purified Chicken anti-Campylobacter Jejuni (206)-Reacts with Fibronectin-binding outer membrane protein.



IF (Immunofluorescence)
Affinity Purified Chicken anti-Campylobacter Jejuni (206)-Reacts with Fibronectin-binding outer membrane protein.

ADDITIONAL INFORMATION

Product Categories/Family for anti-Campylobacter Jejuni antibody	Primary Antibodies; Infectious Agents-Human and Veterinary
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/polyclonal/campylobacter-jejuni-antibody-western-blot-immunohistochemistry-immunofluorescence-elisa/76871>