

GFP/VERO cell line

GFP/VERO Stable Cell Line

CATALOG # AAA78572	APPLICATIONS Flow Cytometry, Functional Assay, Functional Assay
------------------------------	---


PRODUCT SYNONYMS

GFP/VERO; GFP/VERO Stable Cell Line; GFP/VERO cell line

PRODUCT SPECIFICATIONS


Form/Format	Each vial contains 2 ~ 3 x 10 ⁶ cells in 1 ml of 90% FBS + 10% DMSO
Applicable Applications for GFP/VERO cell line	FCM/FACS (Flow Cytometry), FA (Functional Assay)
Preparation and Storage	Immediately upon receipt, store in liquid nitrogen. Ship: Dry ice

IMAGES



FCM/FACS (Flow Cytometry)

Fig-2: Detection of GFP in the GFP/VERO stable cell line through flow cytometry . Parental Vero E6 cells (Green); GFP/VERO cells (Red).



IF (Immunofluorescence)

Fig-1: Analysis of the GFP/VERO stable cell line through fluorescence microscopy. Bright-field image (Left); Fluorescence image (Right).

ADDITIONAL INFORMATION

Related Product Information for GFP/VERO cell line	GFP/VERO Stable Cell Line is a stably transfected Vero E6 cell line which expresses enhanced green fluorescent protein (eGFP). Sequence data: Am eGFPMVSKGEELFTGVVPILVELDGDVNGHKFSVSGEGEGDATYGLTLKFICTTGKLPVPWPTLVTTLTYGVCFSRYPDHMKQHDFFKSAMPEGYVQ
Product Categories/Family for GFP/VERO cell line	Stable Cell Lines

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.