

TLR3/NF-kB LEEPORTER GFP REPORTER CELL LINE | TLR3/NF-kB CELL LINE

TLR3/NF-kB LEEPORTER GFP REPORTER-HEK293 CELL LINE


CATALOG # AAA78590	APPLICATIONS Functional Assay
------------------------------	---

PRODUCT SYNONYMS

TLR3/NF-kB LEEPORTER GFP REPORTER; TLR3/NF-kB LEEPORTER GFP REPORTER-HEK293 CELL LINE; TLR3/NF-kB CELL LINE


PRODUCT SPECIFICATIONS

Form/Format	Each vial contains 2 ~ 3 x 10 ⁶ cells in 1 ml of 90% FBS + 10% DMSO.
Applicable Applications for TLR3/NF-kB cell line	FA (Functional Assay)
Culture conditions	Cells should be grown at 37 degree C with 5% CO ₂ using DMEM medium (w/L-Glutamine, 4.5g/L Glucose and Sodium Pyruvate) supplemented with 10% heat-inactivated FBS and 1% Pen/Strep, plus 3 ug/ml of Puromycin and 15 ug/ml of Blasticidin (Note: Puromycin and Blasticidin can be omitted during the reporter cell assays). It is recommended to quickly thaw the frozen cells upon receipt or from liquid nitrogen in a 37 degree C water-bath, transfer to a tube containing 10 ml of growth medium without Puromycin and Blasticidin, spin down cells, resuspend cells in pre-warmed growth medium without Puromycin and Blasticidin, transfer resuspended cells to T25 flask and culture in 37 degree C-CO ₂ incubator. Leave the T25 flask in the incubator for 1~3 days without disturbing or changing the medium until cells completely recover viability and become adherent. Once cells are over 90% adherent, remove growth medium and passage the cells through trypsinization and centrifugation. At first passage, switch to growth medium containing Puromycin and Blasticidin. Cells should be split before they reach complete confluence. To passage the cells, detach cells from culture vessel with Trypsin/EDTA, add complete growth medium and transfer to a tube, spin down cells, resuspend cells and seed appropriate aliquots of cells suspension into new culture vessels. Subcultivation ration = 1:10 to 1:20 weekly. To achieve satisfactory results, cells should not be passaged over 16 times. Functional validation: A. Response of TLR3/NF-kB LEEPORTER GFP reporter-HEK293 cells to Poly(I:C). 1. Harvest TLR3/NF-kB LEEPORTER GFP reporter-HEK293 cells and seed cells into a tissue culture plate (e.g. 96-well plate in 100 ul of growth medium at 5 x 10 ⁴ cells/well, 24-well plate in 500 ul of growth medium at 2.5 x 10 ⁵ cells/well, 12-well plate in 1 ml of growth medium at 5 x 10 ⁵ cells/well, or 6-well plate in 2 ml of growth medium at 1 x 10 ⁶ cells/well. 2. Incubate cells at 37 degree C in a CO ₂ incubator for overnight. 3. The next day, stimulate cells with different concentrations of Poly(I:C). 4. Incubate at 37 degree C in a CO ₂ incubator for 16 hours. 5. Analyze cells through fluorescence microscopy or flow cytometry.
Preparation and Storage	Immediately upon receipt, store in liquid nitrogen.



Application Data

Fig-2: Induction of TLR3 activity by Poly(I:C) in TLR3/NF-kB LEEporter GFP reporterHEK293 cells was analyzed by flow cytometry. Poly(I:C) was treated at 0.01, 0.05, 0.1, 0.5, 2.5 and 25 ug/ml as noted. Vehicle control (Black line, empty); Poly(I:C) treatment (Green line, grey-filled).



Application Data

Fig-1: Induction of TLR3 activity by Poly(I:C) in TLR3/NF-kB LEEporter GFP reporterHEK293 cells was analyzed by fluorescence microscopy. Poly(I:C) was treated at 0.05, 0.1, 0.5, 2.5 and 25 ug/ml as noted.

ADDITIONAL INFORMATION

Related Product Information for TLR3/NF-kB cell line	The TLR3/NF-kB LEEporter GFP Reporter cell line is a stably transfected HEK293 cell line, which expresses full-length human Toll-like receptor 3 (TLR3) and enhanced green fluorescent protein (eGFP) reporter gene under the transcriptional control of the NF-kB response element. Functional activity of the cell line has been validated by TLR3 ligand assay, in which upon activation by poly (I:C), TLR3 quickly initiates TRIF-dependent signaling pathway and mediates nuclear translocation of NF-kB (Figures 1 and 2).
Product Categories/Family for TLR3/NF-kB cell line	Cell Lines; Reporter Cell Lines

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.