

Human IgG1 Kappa-LE/AF Monoclonal Antibody

Human IgG1 Kappa-LE/AF

CATALOG # AAA78629	HOST Myeloma serum	APPLICATIONS ELISA
------------------------------	------------------------------	------------------------------


PRODUCT SYNONYMS

Human IgG1 Kappa-LE/AF; Monoclonal Antibody; Human IgG1kappa Isotype Control; anti-Human IgG1 Kappa-LE/AF antibody

PRODUCT SPECIFICATIONS

Host	Myeloma serum
Clonality	Monoclonal
Isotype	IgG1 kappa
Purity/Purification	Affinity chromatography on protein G covalently linked to agarose. Purity: >95% as judged by SDS-polyacrylamide gel electrophoresis.
Form/Format	Low Endotoxin, Azide-Free (LE/AF) antibody is supplied as 0.5 mg purified immunoglobulin in 1.0 mL of PBS.
Applicable Applications for anti-Human IgG1 Kappa-LE/AF antibody	FLISA, ELISA
Warning	Reagents contain sodium azide. Sodium azide is very toxic if ingested or inhaled. Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handling. If skin or eye contact occurs, wash with copious amounts of water. If ingested or inhaled, contact a physician immediately. Sodium azide yields toxic hydrazoic acid under acidic conditions. Dilute azide-containing compounds in running water before discarding to avoid accumulation of potentially explosive deposits in lead or copper plumbing.
Preparation and Storage	Aliquot and store at or below -20°C.

IMAGES



ELISA

ELISA plate was coated with serially diluted Human IgG1 Kappa-UNLB. Immunoglobulin was detected with Mouse Anti-Human IgG1 Hinge-BIOT, Mouse Anti-Human IgG2 Fc-BIOT, Mouse Anti-Human IgG3 Hinge-BIOT, Mouse Anti-Human IgG4 pFc'-BIOT, Mouse Anti-Human IgA1-BIOT, Mouse Anti-Human IgA2-BIOT, Mouse Anti-Human IgD-BIOT, Mouse Anti-Human IgE Fc-BIOT, Mouse Anti-Human IgM-BIOT, Mouse Anti-Human Kappa-BIOT, and Mouse Anti-Human Lambda-BIOT followed by Streptavidin-HRP and quantified.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.