

Goat Type I Collagen Antibody | anti-COL1A2 antibody

Goat Anti-Type I Collagen-FITC

CATALOG # AAA78637	HOST Goat	APPLICATIONS ELISA
------------------------------	---------------------	------------------------------


PRODUCT SYNONYMS

Type I Collagen; Goat Anti-Type I Collagen-FITC; anti-COL1A2 antibody

PRODUCT SPECIFICATIONS

Host	Goat
Isotype	IgG
Specificity	Reacts with conformational determinants on type I collagen Referenced species reactivities include – Human; Chicken; Hamster; Mouse; Chinchilla; Porcine; Rat; Elephant; Rabbit ; Bovine 6; Feline; Sheep; Canine; Guinea Pig; Tree Shrew.
Purity/Purification	Affinity chromatography on type I collagen covalently linked to agarose
Sequence Length	397
Applicable Applications for anti-COL1A2 antibody	FLISA, ELISA
Source	Pooled antisera from goats hyperimmunized with type I collagen
Cross Adsorption	Collagen types II, III, IV, V, and VI
Preparation and Storage	Store at 2-8 degree C

IMAGES



Application Data

Human gastric cancer tissue was stained with Goat Anti-Type I Collagen-BIOT followed by Streptavidin-HRP .

References

1. Rasmussen LM, Hansen PR, Ledet T. Homocysteine and the production of collagens, proliferation and apoptosis in human arterial smooth muscle cells. *APMIS*. 2004;112:598-604. (ELISA, Human Reactivity)
2. Kottler UB, Jünemann AG, Aigner T, Zenkel M, Rummelt C, Schlötzer-Schrehardt U. Comparative effects of TGF-beta1 and TGF-beta2 on extracellular matrix production, proliferation, migration, and collagen contraction of human Tenons capsule fibroblasts in pseudoexfoliation and primary open-angle glaucoma. *Exp Eye Res*. 2005;80:121-34. (ELISA, Human Reactivity)
3. Dooley A, Shi-Wen X, Aden N, Tranah T, Desai N, Denton CP, et al. Modulation of collagen type I, fibronectin and dermal fibroblast function and activity, in systemic sclerosis by the antioxidant epigallocatechin-3-gallate. *Rheumatology*. 2010;49:2024-36. (ELISA, WB, Human Reactivity)
4. Zhang WW, Ma KC, Andersen O, Sourander P, Tolleson PO, Olsson Y. The microvascular changes in cases of hereditary multi-infarct disease of the brain. *Acta Neuropathol*. 1994;87:317-24. (IHC-PS, Human Reactivity)
5. Van Deijnen JH, Van Suylichem PT, Wolters GH, Van Schilfgaarde R. Distribution of collagens type I, type III and type V in the pancreas of rat, dog, pig and man. *Cell Tissue Res*. 1994;277:115-21. (IHC-PS, Canine, Human, Porcine, & Rat Reactivity)
6. Moses MA, Wiederschain D, Wu I, Fernandez CA, Ghazizadeh V, Lane WS, et al. Troponin I is present in human cartilage and inhibits angiogenesis. *Proc Natl Acad Sci USA*. 1999;96:2645-50. (IHC-PS, Bovine Reactivity)
7. Dorotka R, Bindreiter U, Macfelda K, Windberger U, Nehrer S. Marrow stimulation and chondrocyte transplantation using a collagen matrix for cartilage repair. *Osteoarthritis Cartilage*. 2005;13:655-64. (IHC-PS, Sheep Reactivity)
8. Weissengruber GE, Egerbacher M, Forstenpointner G. Structure and innervation of the tusk pulp in the African elephant (*Loxodonta africana*). *J Anat*. 2005;206:387-93. (IHC-PS, Elephant Reactivity)
9. Walter I, Schönkypl S. Extracellular matrix components and matrix degrading enzymes in the feline placenta during gestation. *Placenta*. 2006;27:291-306. (IHC-PS, Feline Reactivity)
10. Alcazar O, Cousins SW, Striker GE, Marin-Castano ME. (Pro)renin receptor is expressed in human retinal pigment epithelium and participates in extracellular matrix remodeling. *Exp Eye Res*. 2009;89:638-47. (IHC-PS, WB, Human Reactivity)
11. Chen W, Yu M, Li Y, Chen W, Xia Y. Beneficial effects of astragalus polysaccharides treatment on cardiac chymase activities and cardiomyopathy in diabetic hamsters. *Acta Diabetol*. 2010;47 Suppl 1:S35-46. (IHC-PS, Hamster Reactivity)
12. Kyller M, Witter K, Tichy F. Gingival stippling in dogs: clinical and structural characteristics. *Res Vet Sci*. 2010;88:195-202. (IHC-PS, Canine Reactivity)
13. Zhang Y, Edgley AJ, Cox AJ, Powell AK, Wang B, Kompa AR, et al. FT011, a new anti-fibrotic drug, attenuates fibrosis and chronic heart failure in experimental diabetic cardiomyopathy. *Eur J Heart Fail*. 2012;14:549-62. (IHC-PS, WB, Rat Reactivity)
14. Jiang L, Zhou Y, Xiong M, Fang L, Wen P, Cao H, et al. Sp1 mediates microRNA-29c-regulated type I collagen production in renal tubular epithelial cells. *Exp Cell Res*. 2013;319:2254-65. (IHC-PS, ICC, WB, Rat Reactivity)
15. Sun S, Zhao T, Zhang H, Huang X, Zhang W, Zhang L, et al. Renoprotective effect of berberine on type 2 diabetic nephropathy in rats. *Clin Exp Pharmacol Physiol*. 2015;42:662-70. (IHC-PS, WB, Rat Reactivity)
16. Jerdan JA, Michels RG, Glaser BM. Extracellular matrix of newly forming vessels—an immunohistochemical study. *Microvasc Res*. 1991;42:255-65. (IHC-FS, Guinea Pig & Rabbit Reactivity)
17. Miyaishi O, Sakata K, Matsuyama M, Saga S. Distribution of the collagen binding heat-shock protein in chicken tissues. *J Histochem Cytochem*. 1992;40:1021-9. (IHC-FS, Chicken Reactivity)
18. Dixon IM, Ju H, Reid NL, Scammell-La Fleur T, Werner JP, Jasmin G. Cardiac collagen remodeling in the cardiomyopathic Syrian hamster and the effect of losartan. *J Mol Cell Cardiol*. 1997;29:1837-50. (IHC-FS, Hamster Reactivity)
19. Harada T, Juhn SK, Kim Y, Sakakura Y. Immunohistochemical localization of collagen types I, III, and IV, laminin, fibronectin, and keratin in the endolymphatic sac. *Ann Otol Rhinol Laryngol*. 2000;109:180-6. (IHC-FS, Chinchilla Reactivity)
20. Fraser SA, Crawford A, Frazer A, Dickinson S, Hollander AP, Brook IM, et al. Localization of type VI collagen in tissue-engineered cartilage on polymer scaffolds. *Tissue Eng*. 2006;12:569-77. (IHC-FS, Bovine & Human Reactivity)
21. Chun SY, Lim GJ, Kwon TG, Kwak EK, Kim BW, Atala A, et al. Identification and characterization of bioactive factors in bladder submucosa matrix. *Biomaterials*. 2007;28:4251-6. (IHC-FS, Porcine Reactivity)
22. Albon J, Farrant S, Akhtar S, Young R, Boulton ME, Smith G, et al. Connective tissue structure of the tree shrew optic nerve and associated ageing changes. *Invest Ophthalmol Vis Sci*. 2007;48:2134-44. (IHC-FS, Tree Shrew Reactivity)
23. Behmoaras J, Diaz AG, Venda L, Ko J, Srivastava P, Montoya A, et al. Macrophage epoxygenase determines a profibrotic transcriptome signature. *J Immunol*. 2015;194:4705-16. (IHC-FS, Rat Reactivity)
24. Kothapalli D, Liu S, Bae YH, Monslow J, Xu T, Hawthorne EA, et al. Cardiovascular protection by ApoE and ApoE-HDL linked to suppression of ECM gene expression and arterial stiffening. *Cell Rep*. 2012;2:1259-71. (IHC-FS, ICC, Mouse Reactivity)
25. Lundquist P, Ritchie HH, Moore K, Lundgren T, Linde A. Phosphate and calcium uptake by rat odontoblast-like MRPC-1 cells concomitant with mineralization. *J Bone Miner Res*. 2002;17:1801-13. (ICC, Rat Reactivity)
26. via c-Abl and Egr-1: selective modulation by imatinib mesylate. *Oncogene*. 2009;28:1285-97. (ICC, WB, Human Reactivity)
27. Bhattacharyya S, Ishida W, Wu M, Wilkes M, Mori Y, Hinchcliff M, et al. A non-Smad mechanism of fibroblast activation by transforming growth factor- β . *J Biol Chem*. 2004;279:10333-41. (ICC, WB, Human, Mouse, & Rat Reactivity)
27. Wang J, Chu ES, Chen H, Man K, Go MY, Huang XR, et al. microRNA-29b prevents liver fibrosis by attenuating hepatic stellate cell activation and inducing apoptosis through targeting PI3K/AKT pathway. *Oncotarget*. 2015;6:7325-38. (ICC, WB, Human, Mouse, & Rat Reactivity)

NCBI AND UNIPROT INFORMATION

NCBI GI #	30102
NCBI GeneID	1278
UniProt Accession #	P08123; P02464; Q13897; Q13997; Q13998; Q14038; Q14057; Q15177; Q15947; Q16480; Q16511; Q7Z5S6
NCBI Official Full Name	type I collagen, partial
NCBI Official Synonym Full Names	collagen type I alpha 2
NCBI Official Symbol	COL1A2
NCBI Official Synonym Symbols	OI4
NCBI Protein Information	collagen alpha-2(I) chain
UniProt Protein Name	Collagen alpha-2(I) chain
UniProt Gene Name	COL1A2
UniProt Entry Name	CO1A2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/goat-antibody/type-i-collagen/78637>