

Horse Beta-Endorphin, beta-EP ELISA Kit | beta-EP elisa kit

Horse Beta-Endorphin, beta-EP ELISA Kit

CATALOG #
AAA80174

REACTIVITY
Horse


PRODUCT SYNONYMS

Beta-Endorphin; beta-EP; Horse Beta-Endorphin; beta-EP ELISA Kit; beta-EP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Horse
Specificity	This assay has high sensitivity and excellent specificity for detection of horse beta-EP. No significant cross-reactivity or interference between horse beta-EP and analogues was observed.
Sequence Length	31
Samples	Serum, plasma
Assay Type	Sandwich
Detection Range	3.12 pg/ml-200 pg/ml.
Sensitivity	The minimum detectable dose of horse beta-EP is typically less than 0.78 pg/ml. The sensitivity of this assay, or Lower Limit of Detection (LLD) was defined as the lowest protein concentration that could be differentiated from zero. It was determined the mean O.D value of 20 replicates of the zero standard added by their three standard deviations.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for beta-EP elisa kit	Principle of the assay: This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for beta-EP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any beta-EP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for beta-EP is added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of beta-EP bound in the initial step. The color development is stopped and the intensity of the color is measured.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	119332
UniProt Accession #	P01203
Molecular Weight	3,438 Da
NCBI Official Full Name	Beta-endorphin
UniProt Protein Name	Beta-endorphin
UniProt Gene Name	POMC
UniProt Entry Name	COLI_CAMDR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/horse-elisa-kits/beta-endorphin-beta-ep/80174>