

Mouse leukemia inhibitory factor (cholinergic differentiation factor) ELISA Kit | LIF elisa kit

Mouse Leukemia inhibitory factor, LIF ELISA Kit

CATALOG #
AAA80212

REACTIVITY
Mouse


PRODUCT SYNONYMS

leukemia inhibitory factor (cholinergic differentiation factor); Mouse Leukemia inhibitory factor; LIF ELISA Kit; LIF ELISA kit; CDF; DIA; HILDA; D factor; cholinergic differentiation factor; differentiation inhibitory activity; differentiation stimulating factor; LIF elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of mouse LIF. No significant cross-reactivity or interference between mouse LIF and analogues was observed.
Sequence Length	158
Samples	Serum, plasma, tissue homogenates
Assay Type	Sandwich
Detection Range	15.6 pg/ml -1000 pg/ml.
Sensitivity	The minimum detectable dose of mouse LIF is typically less than 3.9 pg/ml. The sensitivity of this assay, or Lower Limit of Detection (LLD) was defined as the lowest protein concentration that could be differentiated from zero. It was determined the mean O.D value of 20 replicates of the zero standard added by their three standard deviations.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LIF elisa kit	Principle of the assay: This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for LIF has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any LIF present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for LIF is added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of LIF bound in the initial step. The color development is stopped and the intensity of the color is measured.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	88014736
NCBI GeneID	16878
NCBI Accession #	NP_001034626.1
NCBI GenBank Nucleotide #	NM_001039537.2
UniProt Accession #	P09056
Molecular Weight	22,287 Da
NCBI Official Full Name	leukemia inhibitory factor isoform b
NCBI Official Synonym Full Names	leukemia inhibitory factor
NCBI Official Symbol	Lif
NCBI Protein Information	leukemia inhibitory factor; d factor; differentiation-stimulating factor
UniProt Protein Name	Leukemia inhibitory factor
UniProt Gene Name	Lif
UniProt Synonym Gene Names	LIF; D factor
UniProt Entry Name	LIF_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/lif-mouse-elisa-kits/leukemia-inhibitory-factor-cholinergic-differentiation-factor/80212>