

## Human Anti-Ovary Antibody, AOAb ELISA Kit | AOAb elisa kit

Human Anti-Ovary Antibody, AOAb ELISA Kit

CATALOG #  
**AAA80363**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Anti-Ovary Antibody; AOAb; Human Anti-Ovary Antibody; AOAb ELISA Kit; AOAb elisa kit

### PRODUCT SPECIFICATIONS

|                              |                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                   | Human                                                                                                                             |
| Sequence Length              | 96                                                                                                                                |
| Detection Wavelength         | 450 nm                                                                                                                            |
| Sample Volume                | 50-100ul                                                                                                                          |
| Protein Biological Process 1 | others                                                                                                                            |
| Preparation and Storage      | Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions. |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for AOAb elisa kit

**Principle of the Assay:** The microtiter plate provided in this kit has been pre-coated with specificity antigen . Samples are then added to the appropriate microtiter plate wells and incubated. Then add Horseradish Peroxidase (HRP)-conjugated -anti-human IgG to each well and incubate. Finally, substrate solutions are added to each well. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. Calculate the valence of human Anti-Ovary Antibody(IgG) in the samples.

### NCBI AND UNIPROT INFORMATION

|                         |                          |
|-------------------------|--------------------------|
| NCBI GI #               | <a href="#">62363188</a> |
| UniProt Accession #     | <a href="#">Q2KMZ7</a>   |
| Molecular Weight        | 11,324 Da                |
| NCBI Official Full Name | AoaB, partial            |
| UniProt Protein Name    | AoaB                     |
| UniProt Gene Name       | <b>aoaB</b>              |
| UniProt Entry Name      | Q2KMZ7_9NOST             |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/anti-ovary-antibody-aoab/80363>