

Human peroxiredoxin 4 ELISA Kit | PRDX4 elisa kit

Human Peroxiredoxin 4, PRDX4 ELISA Kit

CATALOG #
AAA80721

REACTIVITY
Human


PRODUCT SYNONYMS

peroxiredoxin 4; Human Peroxiredoxin 4; PRDX4 ELISA Kit; Human Peroxiredoxin 4 (PRDX4) ELISA Kit; AOE37-2; PRX-4; thioredoxin peroxidase (antioxidant enzyme) ; PRDX4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Human PRDX4. No significant cross-reactivity or interference between Human PRDX4 and analogues was observed.
Sequence Length	271
Samples	Serum, plasma
Detection Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV% is less than 8% Three samples of known concentration were tested twenty times on one plate to assess. Inter-assay Precision (Precision between assays): CV% is less than 10% Three samples of known concentration were tested in twenty assays to assess.
Detection Wavelength	450 nm
Sample Volume	50-100ul
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PRDX4 elisa kit	Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PRDX4. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific for PRDX4 and Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then a TMB (3,3',5,5' tetramethyl-benzidine) substrate solution is added to each well. Only those wells that contain PRDX4, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The concentration of PRDX4 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	5453549
NCBI GeneID	10549
NCBI Accession #	NP_006397.1
NCBI GenBank Nucleotide #	NM_006406.1
UniProt Accession #	Q13162; Q6FHT3
Molecular Weight	30,540 Da
NCBI Official Full Name	peroxiredoxin-4
NCBI Official Synonym Full Names	peroxiredoxin 4
NCBI Official Symbol	PRDX4
NCBI Official Synonym Symbols	PRX-4; AOE372; AOE37-2; HEL-S-97n
NCBI Protein Information	peroxiredoxin-4; prx-IV; peroxiredoxin IV; antioxidant enzyme AOE372; thioredoxin peroxidase AOE372; thioredoxin peroxidase (antioxidant enzyme); thioredoxin-dependent peroxide reductase AOE372; epididymis secretory sperm binding protein Li 97n
UniProt Protein Name	Peroxiredoxin-4
UniProt Gene Name	PRDX4
UniProt Synonym Gene Names	AOE37-2; Prx-IV
UniProt Entry Name	PRDX4_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/prdx4-human-elisa-kits/peroxiredoxin-4/80721>