

Human progranulin (PGRN) ELISA Kit | PGRN elisa kit

Human progranulin, PGRN ELISA Kit

CATALOG #
AAA80833

REACTIVITY
Human


PRODUCT SYNONYMS

progranulin (PGRN); Human progranulin; PGRN ELISA Kit; Human progranulin (PGRN) ELISA Kit; PGRN elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	104
Samples	Serum, plasma, cell culture supernates
Detection Range	4.7 ng/ml-300 ng/ml
Sensitivity	1.17 ng/ml
Detection Wavelength	450 nm
Sample Volume	50-100ul
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PGRN elisa kit

Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PRGN. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific for PRGN and Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then a TMB (3,3',5,5' tetramethyl-benzidine) substrate solution is added to each well. Only those wells that contain PRGN, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The concentration of PRGN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

268309222

UniProt Accession #	D1LHD5
Molecular Weight	12,259 Da
NCBI Official Full Name	PrgN
UniProt Protein Name	PrgN
UniProt Gene Name	prgN
UniProt Entry Name	D1LHD5_ENTFL

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/progranulin-pgrn/80833>