

Human early prostate cancer antigen 2 (EPCA2) ELISA Kit | EPCA2 elisa kit

Human early prostate cancer antigen 2, EPCA2 ELISA Kit

CATALOG #
AAA80972

REACTIVITY
Human


PRODUCT SYNONYMS

early prostate cancer antigen 2 (EPCA2); Human early prostate cancer antigen 2; EPCA2 ELISA Kit; Human early prostate cancer antigen 2 (EPCA2) ELISA kit; EPCA2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of human EPCA2. No significant cross-reactivity or interference between human EPCA2 and analogues was observed.
Samples	Serum, plasma, tissue homogenates.
Assay Type	Sandwich
Detection Range	12.5ng/ml-400 ng/ml.
Sensitivity	6.25 ng/ml.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<15%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<20%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for EPCA2 elisa kit

Principle of the Assay||This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for EPCA2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells with a Horseradish Peroxidase (HRP) conjugated antibody specific for EPCA2. Following a wash to remove any unbound reagent, a substrate solution is added to the wells and color develops in proportion to the amount of EPCA2 bound in the initial step. The color development is stopped and the intensity of the color is measured.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/early-prostate-cancer-antigen-2-epca2/80972>