

Rabbit Proteasome subunit beta type 9 (PSMB9) ELISA Kit | LMP7/PSMB9 elisa kit

Rabbit Proteasome subunit beta type 9 (PSMB9) ELISA Kit

CATALOG #
AAA82060

REACTIVITY
Rabbit


PRODUCT SYNONYMS

Proteasome subunit beta type 9 (PSMB9); Rabbit Proteasome subunit beta type 9 (PSMB9) ELISA Kit; LMP7/PSMB9 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	50-1000ng/mL
Sensitivity	1.0ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LMP7/PSMB9 elisa kit

Intended Uses: This PSMB9 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Rabbit PSMB9. This ELISA kit for research use only!

Principle of the Assay: PSMB9 ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-PSMB9 antibody and an PSMB9-HRP conjugate. The assay sample and buffer are incubated together with PSMB9-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the PSMB9 concentration since PSMB9 from samples and PSMB9-HRP conjugate compete for the anti-PSMB9 antibody binding site. Since the number of sites is limited, as more sites are occupied by PSMB9 from the sample, fewer sites are left to bind PSMB9-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The PSMB9 concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/lmp7-psmb9-rabbit-elisa-kits/proteasome-subunit-beta-type-9-psmb9/82060>