

Human Intact Parathormone ELISA Kit | I PTH elisa kit

Human Intact Parathormone ELISA Kit

CATALOG #
AAA82686

REACTIVITY
Human


PRODUCT SYNONYMS

Intact Parathormone; Human Intact Parathormone ELISA Kit; I PTH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	1.0-25ng/mL
Sensitivity	0.1ng/mL
Intended Uses	This iPTH ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human iPTH. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for I PTH elisa kit

Principle of the assay: iPTH ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-iPTH antibody and an iPTH-HRP conjugate. The assay sample and buffer are incubated together with iPTH-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the iPTH concentration since iPTH from samples and iPTH-HRP conjugate compete for the anti-iPTH antibody binding site. Since the number of sites is limited, as more sites are occupied by iPTH from the sample, fewer sites are left to bind iPTH-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The iPTH concentration in each sample is interpolated from this standard curve.

Product Categories/Family for I PTH elisa kit
[Signal Transduction](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/i-ptb-human-elisa-kits/intact-parathormone/82686>