

Human Pro-amylin ELISA Kit

Human Pro-amylin ELISA Kit

CATALOG #
AAA83041

REACTIVITY
Human


PRODUCT SYNONYMS

Pro-amylin; Human Pro-amylin ELISA Kit; Pro-amylin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of IAPP. No significant cross-reactivity or interference between IAPP and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between IAPP and all the analogues, therefore, cross reaction may still exist in some cases.
Assay Type	Quantitative Competitive
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Sensitivity	1.0 pg/mL
Intended Uses	This ProlAPP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human ProlAPP. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Pro-amylin elisa kit	Principle of the Assay: IAPP ELISA kit applies the competitive enzyme immunoassay technique utilizing an anti-IAPP antibody and an IAPP-HRP conjugate. The assay sample and buffer are incubated together with IAPP-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the IAPP concentration since IAPP from samples and IAPP-HRP conjugate compete for the anti-IAPP antibody binding site. Since the number of sites is limited, as more sites are occupied by IAPP from the sample, fewer sites are left to bind IAPP-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The IAPP concentration in each sample is interpolated from this standard curve. Intended Uses: This IAPP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human IAPP. This ELISA kit for research use only, not for therapeutic or test applications!
Product Categories/Family for Pro-amylin elisa kit	Signal Transduction

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/pro-amylin/83041>