

## Human Phosphoinositide specific phospholipase C ELISA Kit | PIPLC elisa kit

Human Phosphoinositide specific phospholipase C ELISA Kit

CATALOG #  
**AAA83117**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Phosphoinositide specific phospholipase C; Human Phosphoinositide specific phospholipase C ELISA Kit; PIPLC elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                            |
|-------------------------|----------------------------------------------------------------------------|
| Reactivity              | Human                                                                      |
| Samples                 | Serum, plasma, cell culture supernatants, body fluid and tissue homogenate |
| Sensitivity             | 0.1 ng/mL                                                                  |
| Preparation and Storage | Store all reagents at 2-8 degree C.                                        |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

Related Product Information for PIPLC elisa kit

Intended Uses: This PIPLC ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human PIPLC. This ELISA kit for research use only!

### NCBI AND UNIPROT INFORMATION

|                         |                                                        |
|-------------------------|--------------------------------------------------------|
| NCBI GI #               | 2815230                                                |
| UniProt Accession #     | Q6LAJ0                                                 |
| Molecular Weight        | 30,674 Da                                              |
| NCBI Official Full Name | phosphatidylinositol-specific phospholipase C, partial |
| UniProt Protein Name    | Phosphatidylinositol-specific phospholipase C          |
| UniProt Gene Name       | <b>piplc</b>                                           |
| UniProt Entry Name      | Q6LAJ0_BACTU                                           |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

