

Human Eye muscle antibody ELISA Kit | EM elisa kit

Human Eye muscle antibody ELISA Kit

CATALOG #
AAA83686

REACTIVITY
Human


PRODUCT SYNONYMS

Eye muscle antibody; Human Eye muscle antibody ELISA Kit; EM elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	10-250ng/mL
Sensitivity	1.0ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for EM elisa kit	Intended Uses: This EMAb ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human EMAb. This ELISA kit for research use only! Principle of the Assay: EMAb ELISA kit applies the competitive enzyme immunoassay technique utilizing Eye muscle antigen and an EMAb-HRP conjugate. The assay sample and buffer are incubated together with EMAb-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the EMAb concentration since EMAb from samples and EMAb-HRP conjugate compete for the Eye muscle antigen binding site. Since the number of sites is limited, as more sites are occupied by EMAb from the sample, fewer sites are left to bind EMAb-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The EMAb concentration in each sample is interpolated from this standard curve.
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/em-human-elisa-kits/eye-muscle-antibody/83686>