

Human Corticotropin Like Intermediate Lobe Peptide ELISA Kit | CORT elisa kit

Human Corticotropin Like Intermediate Lobe Peptide ELISA Kit

CATALOG #
AAA83713

REACTIVITY
Human


PRODUCT SYNONYMS

Corticotropin Like Intermediate Lobe Peptide; Human Corticotropin Like Intermediate Lobe Peptide ELISA Kit; CORT elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of CLIP. No significant cross-reactivity or interference between CLIP and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between CLIP and all the analogues, therefore, cross reaction may still exist in some cases.
Assay Type	Quantitative Competitive
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Sensitivity	1.0 pg/mL
Preparation and Storage	Store all reagents at 2-8 degree C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CORT elisa kit	Principle of the Assay: CLIP ELISA kit applies the competitive enzyme immunoassay technique utilizing an anti-CLIP antibody and an CLIP-HRP conjugate. The assay sample and buffer are incubated together with CLIP-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the CLIP concentration since CLIP from samples and CLIP-HRP conjugate compete for the anti-CLIP antibody binding site. Since the number of sites is limited, as more sites are occupied by CLIP from the sample, fewer sites are left to bind CLIP-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The CLIP concentration in each sample is interpolated from this standard curve. Intended Uses: This CLIP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human CLIP. This ELISA kit for research use only, not for therapeutic or test applications!
Product Categories/Family for CORT elisa kit	Human ELISA Kit

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/cort-human-elisa-kits/corticotropin-like-intermediate-lobe-peptide/83713>