

Human Bromodeoxyuridine ELISA Kit | BrDU elisa kit

Human Bromodeoxyuridine ELISA Kit

CATALOG #
AAA84633

REACTIVITY
Human


PRODUCT SYNONYMS

Bromodeoxyuridine; Human Bromodeoxyuridine ELISA Kit; BrDU elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, Plasma, Cell Culture Supernatants, Body Fluid And Tissue Homogenate
Assay Type	Quantitative Competitive
Sensitivity	0.1 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BrDU elisa kit

Intended Uses: This BRD9 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human BRD9. This ELISA kit for research use only!

Principle of the Assay: BRD9 ELISA kit applies the competitive enzyme immunoassay technique utilizing an anti-BRD9 antibody and an BRD9-HRP conjugate. The assay sample and buffer are incubated together with BRD9-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the BRD9 concentration since BRD9 from samples and BRD9-HRP conjugate compete for the anti-BRD9 antibody binding site. Since the number of sites is limited, as more sites are occupied by BRD9 from the sample, fewer sites are left to bind BRD9-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The BRD9 concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/brdu-human-elisa-kits/bromodeoxyuridine/84633>