

Guinea Pig cholesterol ELISA Kit | CH elisa kit

Guinea pig cholesterol ELISA Kit

CATALOG #
AAA84732

REACTIVITY
Guinea Pig


PRODUCT SYNONYMS

cholesterol; Guinea pig cholesterol ELISA Kit; CH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Guinea Pig
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	0.5-10mmol/L
Sensitivity	0.1mmol/L
Intended Uses	This TC ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Guinea pig TC. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CH elisa kit

Principle of the assay: TC ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-TC antibody and an TC-HRP conjugate. The assay sample and buffer are incubated together with TC-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the TC concentration since TC from samples and TC-HRP conjugate compete for the anti-TC antibody binding site. Since the number of sites is limited, as more sites are occupied by TC from the sample, fewer sites are left to bind TC-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The TC concentration in each sample is interpolated from this standard curve.

Product Categories/Family for CH elisa kit

Cardiovascular

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/ch-guinea-pig-pig-elisa-kits/cholesterol/84732>