

## Mouse Heat Shock Protein 70 HSP-70 ELISA Kit | HSP-70 elisa kit

Mouse Heat Shock Protein 70 HSP-70 ELISA Kit

CATALOG #  
**AAA84784**

REACTIVITY  
**Mouse**


### PRODUCT SYNONYMS

Heat Shock Protein 70 HSP-70; Mouse Heat Shock Protein 70 HSP-70 ELISA Kit; Heat Shock Protein 70 HSP-70 (Mouse); HSP-70 elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                    |
|-------------------------|------------------------------------|
| Reactivity              | Mouse                              |
| Assay Type              | Sandwich                           |
| Preparation and Storage | Store all reagents at 2-8 degree C |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for HSP-70 elisa kit

**Intended Uses:** This HSP-70 ELISA kit is intended Laboratory for research use only. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HSP-70 in the sample, this HSP-70 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HSP-70 concentration. The concentration of HSP-70 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

**Principle of the Assay:** This HSP-70 enzyme linked immunosorbent assay applies a technique called a quantitative sandwich immunoassay. The microtiter plate provided in this kit has been pre-coated with a monoclonal antibody specific for HSP-70. Standards or samples are then added to the microtiter plate wells and HSP-70 if present, will bind to the antibody pre-coated wells. In order to quantitatively determine the amount of HSP-70 present in the sample, a standardized preparation of horseradish peroxidase (HRP)-conjugated polyclonal antibody, specific for HSP-70 are added to each well to "sandwich" the HSP-70 immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, A and B substrate solution is added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain HSP-70 and enzyme-conjugated antibody will exhibit a change in colour. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the colour change is measured spectrophotometrically at a wavelength of 450 nm.

Product Categories/Family for HSP-70  
elisa kit

Mouse ELISA Kit

## NCBI AND UNIPROT INFORMATION

|                         |        |
|-------------------------|--------|
| NCBI GI #               | 476003 |
| NCBI Official Full Name | HSP70  |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/heat-shock-protein-70-hsp-70/84784>