

Rat Urokinase ELISA Kit | UK elisa kit

Rat Urokinase ELISA Kit

CATALOG #
AAA85058

REACTIVITY
Rat


PRODUCT SYNONYMS

Urokinase; Rat Urokinase ELISA Kit; UK elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	5.0-100 ng/mL
Sensitivity	1.0ng/ml
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for UK elisa kit

Intended Uses||This UK ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Rat UK. This ELISA kit for research use only!

Principle of the Assay: UK ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-UK antibody and an UK-HRP conjugate. The assay sample and buffer are incubated together with UK-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the UK concentration since UK from samples and UK-HRP conjugate compete for the anti-UK antibody binding site. Since the number of sites is limited, as more sites are occupied by UK from the sample, fewer sites are left to bind UK-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The UK concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

