

Human Anti Teichoic acid antibody ELISA Kit | TA elisa kit

Human Anti Teichoic acid antibody ELISA Kit

CATALOG #
AAA85632

REACTIVITY
Human


PRODUCT SYNONYMS

Anti Teichoic acid antibody; Human Anti Teichoic acid antibody ELISA Kit; TA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	5.0-100ng/mL
Sensitivity	1.0ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TA elisa kit

Intended Uses: This ATAb ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human ATAb.

Principle of the Assay: ATAb ELISA kit applies the competitive enzyme immunoassay technique utilizing Teichoic Acid antigen and an ATAb-HRP conjugate. The assay sample and buffer are incubated together with ATAb-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the ATAb concentration since ATAb from samples and ATAb-HRP conjugate compete for the Teichoic Acid antigen binding site. Since the number of sites is limited, as more sites are occupied by ATAb from the sample, fewer sites are left to bind ATAb-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The ATAb concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

