

## Human 26S proteasome ELISA Kit | 26S PSM elisa kit

Human 26S proteasome ELISA Kit

CATALOG #  
**AAA86468**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

26S proteasome; Human 26S proteasome ELISA Kit; 26S PSM elisa kit

### PRODUCT SPECIFICATIONS

|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| Reactivity                 | Human                                                                      |
| Samples                    | Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate |
| Assay Type                 | Sandwich                                                                   |
| Detection Range            | 5.0-100ng/mL                                                               |
| Sensitivity                | 1.0ng/mL                                                                   |
| Highest Standards          | 100 ng/mL                                                                  |
| Concentration of Standards | Concentration of Standards                                                 |
| Preparation and Storage    | Store all reagents at 2-8 degree C.                                        |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

Product Categories/Family for 26S PSM elisa kit

Cell Biology

### NCBI AND UNIPROT INFORMATION

|                  |           |
|------------------|-----------|
| NCBI GI #        | 284091260 |
| NCBI GeneID      | 8850219   |
| Molecular Weight | 34,833 Da |

|                          |                  |
|--------------------------|------------------|
| NCBI Official Full Name  | 26S proteasome   |
| NCBI Official Symbol     | NAEGRDRAFT_39105 |
| NCBI Protein Information | 26S proteasome   |
| UniProt Protein Name     | 26S proteasome   |
| UniProt Entry Name       | D2VEL9_NAEGR     |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/26s-psm-human-elisa-kits/26s-proteasome/86468>