

Human Heparin associated antibody ELISA Kit | HIT elisa kit

Human Heparin associated antibody ELISA Kit

CATALOG #
AAA87569

REACTIVITY
Human


PRODUCT SYNONYMS

Heparin associated antibody; Human Heparin associated antibody ELISA Kit; HIT elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of HEAb. No significant cross-reactivity or interference between HEAb and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between HEAb and all the analogues, therefore, cross reaction may still exist in some cases.
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Quantitative Competitive or Sandwich
Sensitivity	1.0 ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HIT elisa kit	Intended Uses: This HEAb ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human HEAb. This ELISA kit for research use only, not for therapeutic or test applications! Principle of the Assay: HEAb ELISA kit applies the quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with Heparin Associated antigen. Standards or samples are then added to the microtiter plate wells and HEAb if present, will bind to the Heparin Associated antigen pre-coated wells. In order to quantitatively determine the amount of HEAb present in the sample, a standardized preparation of horseradish peroxidase (HRP)-conjugated HEAb IgG is added to each well to "sandwich" the HEAb immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, substrate solutions are added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain HEAb and enzyme-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The HEAb concentration in each sample is interpolated from this standard curve.
Product Categories/Family for HIT elisa kit	Immunology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/hit-human-elisa-kits/heparin-associated-antibody/87569>