

Mouse Anticardiolipin antibody IgA ELISA Kit | ACA-IgA elisa kit

Mouse Anticardiolipin antibody IgA ELISA Kit

CATALOG #
AAA87853

REACTIVITY
Mouse


PRODUCT SYNONYMS

Anticardiolipin antibody IgA; Mouse Anticardiolipin antibody IgA ELISA Kit; ACA-IgA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Sandwich
Detection Range	5.0-100u/mL
Sensitivity	1.0u/mL
Intended Uses	This ACA-IgA ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Mouse ACA-IgA. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



ADDITIONAL INFORMATION

Related Product Information for ACA-IgA elisa kit

Principle of the assay: ACA-IgA ELISA kit applies the quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with a monoclonal antibody specific for ACA-IgA. Standards or samples are then added to the microtiter plate wells and ACA-IgA if present, will bind to the antibody pre-coated wells. In order to quantitatively determine the amount of ACA-IgA present in the sample, a standardized preparation of horseradish peroxidase (HRP)-conjugated polyclonal antibody, specific for ACA-IgA are added to each well to "sandwich" the ACA-IgA immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, substrate solutions are added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain ACA-IgA and enzyme-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The ACA-IgA concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/aca-iga-mouse-elisa-kits/anticardiolipin-antibody-iga/87853>