

Human Sialic Acid Binding Ig Like Lectin 3 ELISA Kit | SIGLEC3 elisa kit

Human Sialic Acid Binding Ig Like Lectin 3 ELISA Kit

CATALOG #
AAA87899

REACTIVITY
Human


PRODUCT SYNONYMS

Sialic Acid Binding Ig Like Lectin 3; Human Sialic Acid Binding Ig Like Lectin 3 ELISA Kit; SIGLEC3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive or Sandwich
Detection Range	0.5-10ng/mL
Sensitivity	0.1ng/mL
Intended Uses	This SIGLEC3 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human SIGLEC3. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SIGLEC3 elisa kit

Principle of the Assay: SIGLEC3 ELISA kit applies the quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with a monoclonal antibody specific for SIGLEC3. Standards or samples are then added to the microtiter plate wells and SIGLEC3 if present, will bind to the antibody pre-coated wells. In order to quantitatively determine the amount of SIGLEC3 present in the sample, a standardized preparation of horseradish peroxidase (HRP) -conjugated polyclonal antibody, specific for SIGLEC3 are added to each well to "sandwich" the SIGLEC3 immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, substrate solutions are added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain SIGLEC3 and enzyme-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The SIGLEC3 concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/siglec3-human-elisa-kits/sialic-acid-binding-ig-like-lectin-3/87899>