

Monkey 8-Hydroxy-2-desoxyguanosine ELISA Kit | 8-OHdG elisa kit

Monkey 8-Hydroxy-2-desoxyguanosine ELISA Kit

CATALOG #
AAA88273

REACTIVITY
Monkey


PRODUCT SYNONYMS

8-Hydroxy-2-desoxyguanosine; Monkey 8-Hydroxy-2-desoxyguanosine ELISA Kit; 8-OHdG elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Monkey
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Sandwich
Detection Range	5.0-100ng/mL
Sensitivity	1.0ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for 8-OHdG elisa kit

Intended Uses: This 8-OHdG ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Monkey 8-OHdG. This ELISA kit for research use only!

Principle of the Assay: 8-OHdG ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-8-OHdG antibody and an 8-OHdG-HRP conjugate. The assay sample and buffer are incubated together with 8-OHdG-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the 8-OHdG concentration since 8-OHdG from samples and 8-OHdG-HRP conjugate compete for the anti-8-OHdG antibody binding site. Since the number of sites is limited, as more sites are occupied by 8-OHdG from the sample, fewer sites are left to bind 8-OHdG-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The 8-OHdG concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

