

Sheep Nesfatin ELISA Kit | Nesfatin elisa kit

Sheep Nesfatin ELISA Kit

CATALOG #
AAA88642

REACTIVITY
Sheep


PRODUCT SYNONYMS

Nesfatin; Sheep Nesfatin ELISA Kit; Nesfatin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Sheep
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Competitive (Quantitative)
Sensitivity	0.1 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Nesfatin elisa kit

Intended Uses: This NES ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Sheep NES. This ELISA kit for research use only!

Principle of the Assay: NES ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-NES antibody and an NES-HRP conjugate. The assay sample and buffer are incubated together with NES-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the NES concentration since NES from samples and NES-HRP conjugate compete for the anti-NES antibody binding site. Since the number of sites is limited, as more sites are occupied by NES from the sample, fewer sites are left to bind NES-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The NES concentration in each sample is interpolated from this standard curve.

Product Categories/Family for Nesfatin elisa kit

Neurobiology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

