

Rat Hyperglycosylated chorionic gonadotropin ELISA Kit | hCG-H elisa kit

Rat Hyperglycosylated chorionic gonadotropin ELISA Kit

CATALOG #
AAA89772

REACTIVITY
Rat


PRODUCT SYNONYMS

Hyperglycosylated chorionic gonadotropin; Rat Hyperglycosylated chorionic gonadotropin ELISA Kit; hCG-H elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Sandwich
Detection Range	50-1000nmol/L
Sensitivity	1.0nmol/l
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for hCG-H elisa kit	For samples: Serum, plasma, cell culture supernatants, body fluid and tissue homogenate INTENDED USE This hCG-H ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Rat hCG-H. This ELISA kit for research use only! PRINCIPLE OF THE ASSAY hCG-H ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-hCG-H antibody and an hCG-H-HRP conjugate. The assay sample and buffer are incubated together with hCG-H-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the hCG-H concentration since hCG-H from samples and hCG-H-HRP conjugate compete for the anti-hCG-H antibody binding site. Since the number of sites is limited, as more sites are occupied by hCG-H from the sample, fewer sites are left to bind hCG-H-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The hCG-H concentration in each sample is interpolated from this standard curve.
---	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/hcg-h-rat-elisa-kits/hyperglycosylated-chorionic-gonadotropin/89772>