

## Chicken Glial Fibrillary Acidic Protein ELISA Kit | GFAP elisa kit

Chicken Glial Fibrillary Acidic Protein ELISA Kit

CATALOG #  
**AAA89840**

REACTIVITY  
**Chicken**


### PRODUCT SYNONYMS

Glial Fibrillary Acidic Protein; Chicken Glial Fibrillary Acidic Protein ELISA Kit; GFAP elisa kit

### PRODUCT SPECIFICATIONS

|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| Reactivity                 | Chicken                                                                    |
| Assay Type                 | Competitive or Sandwich                                                    |
| Samples                    | Serum, plasma, cell culture supernatants, body fluid and tissue homogenate |
| Sensitivity                | 0.1 ng/mL.                                                                 |
| Highest Standards          | 50ng/mL                                                                    |
| Concentration of Standards | Concentration of Standards                                                 |
| Preparation and Storage    | Store all reagents at 2-8 degree C                                         |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for GFAP elisa kit | <p><b>For samples:</b> Serum, plasma, cell culture supernatants, body fluid and tissue homogenate</p> <p><b>INTENDED USE</b> This GFAP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Chicken GFAP. This ELISA kit for research use only!</p> <p><b>PRINCIPLE OF THE ASSAY</b> GFAP ELISA kit applies the quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with a monoclonal antibody specific for GFAP. Standards or samples are then added to the microtiter plate wells and GFAP if present, will bind to the antibody pre-coated wells. In order to quantitatively determine the amount of GFAP present in the sample, a standardized preparation of horseradish peroxidase (HRP)-conjugated polyclonal antibody, specific for GFAP are added to each well to “sandwich” the GFAP immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, substrate solutions are added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain GFAP and enzyme-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The GFAP concentration in each sample is interpolated from this standard curve.</p> |
| Product Categories/Family for GFAP elisa kit   | Neurobiology                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

NCBI AND UNIPROT INFORMATION

|                                  |                                                                                        |
|----------------------------------|----------------------------------------------------------------------------------------|
| NCBI GI #                        | 251802                                                                                 |
| NCBI GeneID                      | 2670                                                                                   |
| UniProt Accession #              | P14136; Q53H98; Q5D055; Q6ZQS3; Q7Z5J6; Q7Z5J7; Q96KS4; Q96P18; Q9UFD0; B2RD44; D3DX59 |
| Molecular Weight                 | 49,880 Da                                                                              |
| NCBI Official Full Name          | glial fibrillary acidic protein                                                        |
| NCBI Official Synonym Full Names | glial fibrillary acidic protein                                                        |
| NCBI Official Symbol             | GFAP                                                                                   |
| NCBI Protein Information         | glial fibrillary acidic protein                                                        |
| UniProt Protein Name             | Glial fibrillary acidic protein                                                        |
| UniProt Gene Name                | GFAP                                                                                   |
| UniProt Synonym Gene Names       | GFAP                                                                                   |
| UniProt Entry Name               | GFAP_HUMAN                                                                             |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/gfap-chicken-elisa-kits/glial-fibrillary-acidic-protein/89840>